




www.venturemedical.com

Venture Medical ReQuip, Inc.  
6008 Bonacker Drive  
Tampa, FL 33610  
800-627-3215 | 813-740-8640  
Fax: 800-466-1251  
sales@venturemedical.com

## Stryker 1089 Eye Surgery Stretcher - Refurbished



- New Patient Pad
- 8" wheels
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual-side hydraulic controls
- Fold-down side rails
- Four IV receptacles
- Central-locking brakes
- Foot-end brake and steer pedal
- Hands-free Trendelenburg
- Integrated storage tray with oxygen bottle holder
- Retractable fifth-wheel steering
- Rubber bumpers
- Stationary foot-end litter
- Uni-lower pedal
- Dual-articulation, enhanced clearance headpiece

The Stryker 1089 Eye Surgery Stretcher offers maximum head-end clearance, giving surgeons ample leg room to maneuver while performing delicate eye procedures. The reliable work surface and mobile design combine to help with surgical throughput.

The Stryker 1089 features Positioning controls that tuck neatly within the headpiece and allow full articulation. Enhanced clearance and dual-articulating headpieces maximize surgical access.

The Stryker 1089s Flexible tube doubles as air delivery and drape support: Dual cylinder pneumatic backrest, Dual pedestal hydraulics, and Dual side hydraulic controls.

### Specifications:

- Model Number: 1089
- Overall Length: 90"
- Overall Width (side rails up): 31.5"
- Weight Capacity: 500 lbs.
- Height Range (to litter top): 22.25" to 34"
- Littler Positioning:
  - Backrest: 0° to 90°
  - Knee Gatch: 0° to 30°
- Trendelenburg/Reverse Trendelenburg.:  $\pm 18^\circ$
- Patient Surface: 26" x 87"
- Siderails: 13" x 55"
- Caster Diameter: 8

**VENTUREMEDICAL.COM**

© 2007-2025 Venture Medical ReQuip, Inc.