

Getinge 5624 Warming Cabinet Refurbished



Stock Photo – Options May Vary

The Getinge 5624 warming cabinet is a great choice for warming blankets and fluids when space is at a premium, but a double door unit is required.

The 5624 features glass door panes which enables the staff to view what is inside the unit without having to open the door and lose heat in the process. Digital displays and controls make this modern design very user friendly.

The Getinge 5624 warming cabinet comes completely refurbished with a comprehensive warranty.

- Full-Size, 2-Door Cabinet
- Compact Design
- Digital Display and Controls
- Standard 120V Plug
- Stainless Steel Construction
- Glass Door Panes
- Over Temp Alarms
- Completely Refurbished

Specifications:

- Exterior Dimensions:
30"W x 60"H x 24.7"D
- Interior Dimensions (Upper):
26"W x 13"H x 22.7"D
- Interior Dimensions (Lower):
26"W x 40"H x 22.7"D
- Capacity: 4.1 Cubic Feet Upper
- Capacity: 13.4 Cubic Feet Lower
- Power: 115V, 60 Hz
- Current Draw: 12 Amps
- Temp Range: 90° - 160°
- Item Weight: 303 Lbs

***Brochure and images are a general overview of the product.
Actual product options may vary.*