

# power200

power exam table

## Power Hi-Low, Manual Back Exam Table with OneTouch Patient Positioning® System, Hand and Foot Control

### For examinations

The power200 is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers.

Designed for exams, the power200 meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.



Model 4070-650-200

**PROTECTS** > patients from infection > facilities from damage > staff from injury

### Key features at a glance

- > Hand Control with OneTouch Patient Positioning®, 2 programmable positions & return-to-home
- > Power Hi-Lo with smooth moving articulating back adjustment
- > One function foot control
- > Dual electro-mechanical actuator system for efficient speed
- > 18"-21" paper roll holder (53.3cm)
- > 650lb (294.84kg) patient weight capacity
- > Low-access ADA height 18.5" (±0.5") - 46.9cm (±1.27cm)
- > Pelvic tilt (7°)
- > Drain Pan: stainless steel 1.5" deep (3.81cm)
- > Front storage drawer
- > Side storage drawer – reversible
- > Removable, seamless upholstery with Prefixx® Protective Finish
- > Stirrups with four lateral positions, adjustable length
- > Front leg cushion for patient comfort
- > 3 year limited warranty
- > Heavy gauge steel extension leg pull-out to accommodate all patient heights
- > Hospital grade receptacle 110V with child safety cover (patient left)
- > Durable galvanized steel construction
- > Easy-clean powder coat enamel surfaces

## Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm) at seat max. 37" at seat (93.9cm)
- › Width: 28" (71.12cm)
- › Length: 57" – 76" with foot section extended (144.7cm - 193cm)
- › Leg Rest: 19" W x 13" D (48.26cm W x 33cm D)
- › Front storage drawer: 19" W x 12" D (48.26cm W x 30.48cm D)
- › Stirrups: four lateral positions, adjustable length
- › Dual electro-mechanical actuator system for efficient speed
- › Electrical: 120 VAC, 60 HZ
- › International model: 4070-650-201 220/240V, 50/60Hz
- › Net weight: 440lb (200kg)
- › Shipping weight: 480lb (218kg)
- › Shipping dimensions: L64 x W30 x H31 (163cm x 76cm x 79cm)

## Accessories

- › 042 – Patient assist arms complete (pair) field installable
- › 045 – Universal Accessory Rail-Seat Section
- › 234 – Two Hook IV pole kit (includes IV pole and clamp)
- › 235 – Four Hook IV pole kit (includes IV pole and clamp)
- › 236 – Fixed armboard package (EA)
- › 237 – Articulating armboard package (EA)
- › 251 – Bierhoff knee crutches (pair)
- › 467 – Side Mounting plate for Welch Allyn Light Bracket #48955 (Patient Right), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 468 – Side Mounting plate for Welch Allyn Light Bracket #48955 (PatientLeft), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 471 – Front Mounting Plate for Welch Allyn Light Bracket #48955 (Patient Left), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 472 – Front Mounting Plate for Welch Allyn Light Bracket #48955 (Patient Right), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 822 – Removable stainless steel drain pan W/spout
- › 870 – 110v drawer warmer (factory installed)
- › 872 – 220v drawer warmer (factory installed) (international)
- › 910 – Easy power table installation kit

## Colors



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

©2020 United Metal Fabricators Inc. dba UMF Medical 05.23 All Rights Reserved. OneTouch Patient Positioning® is a registered trademark of United Metal Fabricators Inc. 5-1335LT Rev A