

### Infection Control • Patient Comfort • Efficiency



### Color Chart and Part Number:

<b>4070-650RI</b> River Rock	<b>4070-650SN</b> Sand Grey	<b>4070-650MO</b> Mocha	<b>4070-650AD</b> Adobe
<b>4070-650ST</b> Steel Blue	<b>4070-650MI</b> Midnight Blue	<b>4070-650SG</b> Sage	<b>4070-650ON</b> Onyx

### Standard Features

- Power Hi-Lo with air-spring assisted manual back adjustment
- Low access ADA compliant 18.5" (± .05") (46.9cm)
- OneTouch WheelBase™ for infection control
- 650 lb (294.84 kg) patient weight capacity
- Magnetic hand control with OneTouch Patient Positioning™—two programmable table positioning buttons, plus return-to-home
- Two function foot control with hi-low functions
- Large front storage drawer and reversible rear storage drawer
- Dual electro-mechanical actuator system for efficient speed
- 12 gage steel extension leg-pull out to accommodate for all patient heights
- Removable upholstered top with seamless antimicrobial guard
- Seamless upholstered leg cushion elevates 10° for patient comfort
- Paper roll holder (18" - 21") (45.7cm-53.3cm)
- One touch, easy glide stirrups with 4 lateral positions and friction lock length adjustments
- Stainless steel drain pan (1-1/2" deep) (3.81cm)
- Hospital grade receptacle 110V with child safety cover (patient left)
- Drawer warmer 110V
- Pelvic tilt (7°)
- 3 year limited warranty

### Optional Accessories

- **014** - Leg section booster pad (field installable)
- **251** - Bierhoff knee crutches (pair)
- **822** - Stainless steel drain pan with spout
- **823** - Deep stainless steel 4" (10.1cm) drain pan with spout
- **848** - Drawer warmer 220V (factory installed) with child safety cover (patient left)
- **856** - 220/240 VAC, 50/60HZ (factory installed)
- **4040-650-300F** - Fire retardant foam (CAL133)

### Table Specifications

- **Height:** Min. 18.5" at seat (± .05") (46.9cm)  
Max. 37" at seat (93.9cm)
- **Width:** 27" (68.6cm)
- **Length:** 57" - 76" with foot section extended (144.7cm - 193cm)
- **Leg rest:** 21"W x 13"D (53.3cm W x 33cm D)
- **Front Storage Drawer:** 19"W x 12"D
- **Rear Storage Drawer:** 12.5"W x 24.5"D
- **Electrical:** 120 VAC, 60 HZ or 220/240VAC, 60HZ
- **Stirrups:** Four lateral positions/infinately adjustable