


OptiVac® G180

AC/DC Portable Suction Unit



G180 shown with Disposable Hydrophobic/Bacteria Filter and optional Carrying Case

The Quality Of Gomco® Where You Need It Most.

The G180 is the second generation OptiVac® in the industry-leading Gomco line of portable suction equipment. Its compact design, ease of operation, rechargeable battery, AC power supply, and fully-enclosed unit are just a few of its advanced features.

FEATURES

- Run time operation on fully-charged battery for approximately 3 hours at maximum vacuum.
- 25-550 mm Hg vacuum range, coupled with greater than 30 lpm flow, provides powerful suction for faster aspiration.
- Fully chargeable in 8 hours.
- Will run on AC power even with a low battery.
- Will run for approximately 1 hour after being charged for just 2 hours.
- Long-life rechargeable battery eliminates the need for frequent replacement.
- Sealed water-resistant unit is easy to clean, and minimizes potential contamination or damage to the compressor and battery.
- Control panel at an angle for easy viewing.
- Quality-engineered thermoplastic provides impact resistance and durability.
- 3-year warranty on the unit, and 1 year on the battery.





Disposable Hydrophobic/Bacteria Filter



Suction Tubing Kit



Disposable Collection Canister

SPECIFICATIONS

Weight:	11.4 lbs. with internal power
Dimensions:	9.4" H x 7.5" W x 16.8" L
Vacuum Range:	25-550 mm Hg
Vacuum Flow:	>30 lpm
Battery:	Rechargeable lead acid
Battery Life (at maximum vacuum):	3 hours at 550 mm Hg
Battery Charge Time:	Full charge in 8 hours
Battery Depleted (unit automatically shuts off):	Yes
Bad Battery Sensor:	5 red LEDs
Battery Charge Status:	Indicator lights
Remote/Internal AC Charger:	Internal
Charging Circuit Capable of Running on Low Battery:	Yes
Collection Canister:	1500 ml
AC Operation:	110-230 volts, 50/60 Hz
DC Operation:	12 volt DC battery

Classified per UL 60601-1 and U.L.C.

ORDERING INFORMATION

G180	OptiVac® Portable Suction Unit with Gauge & Regulator, ON & OFF Buttons
L190-BAG	Carrying Case with shoulder strap
20-08-0003	Disposable Collection Canister, 1500 ml, stem inlet, 48/case
20-08-0016	Disposable Collection Canister, 1500 ml, stem inlet, 16/case
01-90-3928	Disposable Hydrophobic/Bacteria Filter, 1/4" x 1/4" barb, 3 filters/pkg.
S610100	Suction Tubing Kit (1 each 72" length and 15" length, 1/4" ID, clear, PVC)

Tubing can be cut to desired length and attached to disposable hydrophobic/bacteria filter.



Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, Missouri 63110
 314/771-2400 • 800/444-3954 • Fax: 314/771-0650 • www.alliedhpi.com
 International offices in Europe, Asia, and South America

U.S.

Toll Free: 800/444-3940
 Fax: 800/477-7701

©Allied Healthcare Products, Inc.

Canada & Mexico

Toll Free: 800/446-0552
 Fax: 800/246-6201

International

Telephone: 314/268-1683
 Fax: 314/771-5183

Form No. 21-00-0002 12/04