



STRETCHER

## Designed for Outstanding Mobility & Safety

A feature-rich standard configuration to meet your needs and your budget.



**Fold-away push handles**  
improve maneuverability



**Sliding Backrest**  
reduces patient migration resulting in less boosting



**Head of bed angle indicator**  
aids with airway clearance head elevation compliance.



**Knee gatch**  
provides the patient added comfort. Easy to use manual crank.



**Fifth wheel**  
allows optimal handling and control during transport.

**Three-year  
Limited Product  
Warranty**

**700 lbs Safe  
Working Load**

# STRETCHER



## Examination and Recovery Configuration

Includes Knee Gatch  
Recovery Position

## Transport Configuration

Stationary foot-end  
patient platform

### Standard Features

- 700 lbs weigh capacity
- 28-inch width patient deck
- Fifth wheel
- Knee Gatch /Recovery Position (on Examination and Recover configuration only \*)
- Fold-away push handles on head-end
- Sliding backrest
- Permanently mounted IV pole at head-end, 2 hooks
- Angle indicators for head of bed and trend
- 8-inch integral casters with wheel covers
- Pneumatic assisted backrest
- Trendelenburg and reverse trendelenburg
- Central-locking brake system
- Brake and steer pedals on head and foot ends
- Hydraulic elevation control pedals on both sides
- Fold down siderails
- 4 integrated IV sockets
- Drainage bag holders
- Bumper wheels on each corner
- Patient storage tray on base
- Oxygen tank holder on base
- 4" Comfort Series Pad

### Specifications

- Weight Capacity 700 lb
- Patient surface width 28"
- Patient surface length 77"
- Knee gatch 0° - 45
- Backrest 0° - 90°
- Height range 23.25" - 34"
- Trendelenburg 17°
- Reverse trendelenburg -17°
- Caster diameter 8"
- Overall length 84"
- Overall width 36"

### Optional Features

- Deck Mounted Oxygen Tank Holder
- Fold-away push handles on foot-end
- Side Rail Pads
- Chart holder
- Removable IV Pole, 2 stage and 2 hooks

### Warranty

Three-year limited product warranty

MCM9010	Premium Transport/Recovery Stretcher, Knee Gatch, IV Pole, Base Oxygen Tank Holder, 700 lbs. weight capacity
MCM9022	Foldable Push Handles on Foot End
MCM9023	Removable IV Pole, adjustable height, 2 Hooks
MCM9024	Chart Holder
MCM9025	Side Rail Pads
MCM9026	Deck Mounted Oxygen Tank Holder