

## Medco Aspirator IV Liposuction Unit – New



*Stock Photo – Options May Vary*

The Aspirator IV utilizes two advanced, state of the art positive displacement pumps. The aspirators' dual headed pumps relieve back pressure to allow each individual pump to pull against the patient so that a deep vacuum is quickly reached and maintained. If one pump should go down, the other three will continue to function independently.

The positive displacement pumps used in Medco aspirators feature the quietest on the market. The aspirator is made even quieter with 1" shock mounts to eliminate vibration. This, combined with a super quiet muffler and insulation system, make the Medco Aspirator IV the quietest aspirator on the market today.

The Medco Aspirator IV is AC-powered, while the foot pedal is pneumatic and not powered by electricity, per OSHA standards

The Medco Aspirator comes with a 5 Year Mfg. Warranty.

- Made in the USA
- AC- Powered
- Quietest on the Market
- Insulation System
- Dual Headed Pump
- Pneumatic Foot Pedal
- 5 Year Warranty

### Specifications:

- Overall Dims :
  - 17"W x 18"H x 18"D
- Weight: 64 lbs
- Voltage: 110 or 220 V.
- Max Vacuum: 29+” Hg.
- Max Flow: 7.5 CFM (212.38LPM)
- Sound Level: 29± 3 db.
- Pump: 2 x 2 Piston
- Horsepower: 2/3 hp.
- Warranty: 5 Year Mfg.

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*