

# Passport® 8

## Specifications

### Physical Dimensions

Monitor Size: 22.4 cm (H) x 13.5 cm (D) x 24.2 cm (W)  
8.8" (H) x 5.3" (D) x 9.5" (W)

Weight: 3.5 kg (7.8 lbs) (Standard parameters configuration, including a lithium battery and a recorder)

### Input/Output Communications ECG Analog Output Specification

Bandwidth (-3dB; reference frequency: 10Hz)  
Diagnostic Mode: 0.05 to 150 Hz  
Monitor Mode: 0.5 to 40 Hz  
Surgical Mode: 1 to 20 Hz  
Monitor QRS Delay: ≤ 25 ms (in diagnostic mode, and non-paced)  
Sensitivity: 1V/mV ±5%

### Pacer Enhancement

Amplitude: Voh ≥ 2.5V  
Width: 10 ms ± 5%  
Signal Rising and Falling Time: ≤ 100 μs

### Communication Protocol

DIAP Communication Protocol  
HL7 Communication Protocol

### Electrical Ratings

Power Supply: AC power, battery

### AC Mains Power Source

AC Power: 100 to 240 VAC (±10%)  
Current: 1.3 to 0.5 A  
AC Power Input Frequency: 50/60 Hz (±3Hz)

### Environmental Conditions

#### Temperature Range

Operating: 10°C to 40°C  
Non-operating (storage): -20°C to 60°C

#### Humidity Range

Operating: 15% to 95%, non-condensing  
Non-operating (storage): 10% to 95%, non-condensing

#### Altitude Range

Operating: (-1,000 to 9,841 feet ASL)  
1,074 hPa to 700 hPa  
(806 mmHg to 394 mmHg)  
Non-operating (storage): (-1,000 to 15,000 feet ASL)  
1,074 hPa to 700 hPa



(806 mmHg to 394 mmHg)

### Monitor Interface

1 AC Power Input  
Connectors:  
1 RJ45 connector, 100 Base-TX, IEEE 802.3  
1 VGA connector  
2 USB 2.0  
1 Equipotential Grounding Terminal  
1 Multifunctional connector: Outputs defibrillator synchronization signals, nurse call signals and analogy output signals


### Battery

Type: Chargeable Lithium-Ion, 11.1DVC, 4.5 Ah  
Number of Batteries: 1  
Run Time: 6 hours when powered by a new fully-charged battery (25°C, SpO2 sensor and ECG cable connected, Temp cable not connected, Auto NIBP measurements at an interval of 15 minutes)  
Charge Time: Less than 3 hours to 90%, and less than 4 hours to 100% when the monitor is off.  
Less than 8 hours to 90%, and less than 12 hours to 100% when the monitor is on.  
Shutdown Delay: at least 20 min (after a low battery alarm first occurs)

Top

Front

Side




# Passport® 8

## Specifications

### Mounting Solutions


#### Variable Height Mount (VHM)

Moment Load: 40 lbs (18.1 kg)



#### Pivot Arm Wall Mount

Moment Load: 60 lbs (27.2 kg)




### Rolling Stand

Dimension: 47" (H) x 27" (W) x 27" (D)  
119.4 cm (H) x 68.6 cm (W)  
x 68.6 cm (D)

Moment Load: 35 lbs (15.9 kg)

Max Head Tilt: 10 degrees



Mindray DS USA, Inc.  
800 MacArthur Blvd., Mahwah, NJ 07430  
Tel: 800.288.2121 Tel: 201.995.8000 Fax: 800.926.4275 [www.mindray.com](http://www.mindray.com)

©2015 Mindray DS USA, Inc.  
Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd.  
03/15 P/N: 0002-08-8002 Rev B

**mindray**  
NORTH AMERICA