

Hill-Rom Advanta P-1600 Hospital Bed - Refurbished



Stock Photo – Options May Vary

The Hill-Rom Advanta is one of the most popular full electric, adjustable hospital beds. The Advanta is great for home use, nursing, schools, and medical facilities.

The Hill-Rom Advanta features promote safety and simplicity. It includes Easy Grip siderail controls that are conveniently located in a central location where patient and caregiver interaction occurs.

The Hill-Rom Advanta comes completely refurbished with a comprehensive warranty.

- Head Section Inclination
- Automatic Contour
- Articulating Frame
- Knee Section Inclination
- Foot Section Inclination
- Trendelenburg/Rev Trend
- Chair Position
- CPR Release
- Central Brake and Steer
- Manual Operation
- Headboard
- Footboard
- Completely Refurbished

Specifications:

- Overall Length: 91" (231 cm)
- Retracted Length: 78.5" (199.4 cm)
- Mattress Support Length: 83" (211 cm)
- Maximum Width (siderails up): 42.75" (108.59 cm)
- Maximum Width (siderails stored): 36" (91 cm)
- Maximum Headboard Height: 48" (122 cm)
- Maximum Footboard Height: 45" (114 cm)
- Maximum Siderail Height: 45" (114 cm)
- Head Section Inclination (maximum): 65°
- Knee Section Inclination (maximum): 35°
- Automatic Contour Up (maximum): 15°
- Automatic Contour Down: 28°
- Automatic Contour: 10°
- Hi-low Range: 13.75" (33.66 cm)
- Trendelenburg: 11°
- Reverse Trendelenburg: 11°
- Maximum Patient Weight: 500 lb (226.8 kg)

***Brochure and images are a general overview of the product. Actual product options may vary.*