

Hill-Rom VersaCare Hospital Bed - Refurbished



Stock Photo – Options May Vary

The Hill-Rom VersaCare med-surg bed is a state-of-the-art bed that integrates the high quality of Hill-Rom products with sophisticated technology advancements for both the patient and the caregiver. Its versatility is unmatched.

The Point-of-Care control panel is easy to read and use, using simple icons for bed adjustment and control. The *HandsFree* foot pedal system allows the caregiver to adjust bed height and the backrest without using their hands.

The VersaCare is designed to be adjusted to a very low height, which assists patients getting in and out of the bed. Safeview alerts also show a visual status to indicate current conditions of the bed, including status of side rails, brakes, bed height, and more. If the head of the bed is lowered more than 30°, the Head-of-Bed alarm will alert the caregiver that attention is needed. Optional integrated bed exit alarms detect whether the patient is in bed, sitting on the edge or has exited.

The Hill-Rom Versacare comes completely refurbished with a mattress and comprehensive warranty.

- New Foam Mattress
- Head/Knee/Hi-Lo Adjustable
- Trend/Rev Trend
- “Brake Not Set” Alarm
- FlexAFoot Mechanism
- Point-Of-Care Controls
- Maintains Low Height in Chair Position
- Head/Foot Boards
- Chair Position
- CPR Release
- Central Brake and Steer

Specifications:

- Length: 94.5" (2400 mm)
- Width: 40" (1020 mm)
- Minimum Bed Height: 18" (457 mm)
- Maximum Bed Height: 37" (940 mm)
- Mattress Length: 86" (2180 mm)
- Mattress Width: 35.5" (902 mm)
- Maximum Head Elevation: 0° - 65°
- Maximum Knee Elevation: 0° - 16°
- Maximum Trend/Rev Trend: 15° / 10°
- Maximum Patient Capacity: 500 lbs. (227 kg)
- Safe Working Load: 550 lbs
- Caster Size: 5" (127 mm)

***Brochure and images are a general overview of the product. Actual product options may vary.*