

Mindray Passport 12m Patient Monitor Refurbished



Stock Photo – Options May Vary

The Mindray Passport 12m was designed with simplicity and user-friendliness in mind. The large (12.1”) touchscreen display has auto adjustable numerics and can display up to 8 waveforms. The user-friendly design and menus make the system easy to use and understand.

Included with the unit: ECG Cable (3 Lead), Reusable SpO2 finger sensor, NIBP Hose and Adult BP Cuff, Roll of Recorder Paper, New Batteries, Power Cord, and Users Manual. The monitor comes fully refurbished with a comprehensive warranty.

- Large, Bright 12.1” Color Display
- Touchscreen
- Measures ECG, NIBP, SpO2, Temp
- 3 or 5 Lead ECG
- SpO2
- Onboard Thermal Printer
- Fully User Configurable Display
- 3 to 8 Waveforms Displayed
- Trending
- EMR Compatible
- Battery Operated w/ Built-In Charger
- Lightweight Design
- Carry Handle
- Mounting Options Available
- Full Range of Accessories Available
- Optional Co2 Available
- Optional Anesthesia Gas Available
- Completely Refurbished

Specifications:

- Monitor Dims: ~11.7”H x 13.2” W x 7.4” D
- Display Size (Screen): 12.1”
- Display Type: TFT LCD Color
- Recorder Paper: 50mm, Thermal
- Battery: Rechargeable, Lith Ion
- Battery Run Time: Approx 300 minutes
- Weight: Approx 16.5 lbs
- Power Requirement: 120v

***Brochure and images are a general overview of the product. Actual product options may vary.*