

proglide4010
head centric procedure chair



The Next Generation Head Centric Procedure Chair just got better



Shown with optional Armrests

The ONLY Power Procedure Chair with Wheels and ADA Access
Helps you protect patients from infection
Helps give you maximum room flexibility
Helps your staff to be more efficient



CHAIR SPECIFICATIONS

Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at seat Max. 37" at seat (93.9cm)
Seat Width: 28" (71.12cm)
Back Width: 24.5" (62.23cm)
Length: 70" – 80" with foot section extended (177.8cm – 203.2cm)
Leg Rest: 28" W x 13" D (71.12cm W x 33.02cm D)
Dual Electro-Mechanical Actuator System for Efficient Speed
Electrical: 120 VAC, 60HZ
International Model 4010-650-301: 220/240V, 50/60Hz

ACCESSORIES

047 Adjustable/self leveling arms pair- field installable
 050 Removable U-shaped headrest
 051 Removable rectangular headrest
 467 Mounting plate for Welch Allyn Light #48955 - patient right
 468 Mounting plate for Welch Allyn Light #48955 - patient left
 830 Patient safety straps
 900 Lifting bracket kit
 920 Easy installation kit
 Fire retardant foam (CAL133)) available

KEY FEATURES

model # 4010-650-300

Chair Designed for Head and Neck Centric Practices
 OneTouch Wheelbase®
 650 lb. (294.84 kg) Patient Weight Capacity
 Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
 4 Power Functions: Hi-Lo, Back, Tilt, Leg Lift & Extension
 Hand Control with OneTouch Patient Positioning®
 2 Programmable Positions & Return-to-Home
 Foot Control with OneTouch Patient Positioning®
 2 Programmable Positions and Return-to-Home
 Double-articulating Headrest with Height Adjustability
 and Interchangeable Headrest Options
 Slim Back Cushion for Accessibility
 Removable Seamless Upholstered Cushion
 Trendelenburg Position
 Hospital-Grade Electrical Receptacle with Child Safety Cover
 (Patient Left)
 Concealed Paper Roll Holder
 3 Year Limited Warranty

proglide4010

head centric procedure chair



COLORS

Shown with optional Armrests (047)



Why wouldn't you want a procedure chair that moves?

Features that matter - Empowering the Patient Experience

Infection Control

One Touch WheelBase® for cleaning surfaces around and under the Procedure Chair. Seamless upholstery and smooth body surfaces make it fast and easy to wipe clean. No cracks or crevices for mold or bacteria to grow. The upholstered top is removable for deep cleaning and replacement to extend the life of the procedure chair.

OneTouch Patient Positioning®

Personalize comfort and ergonomics with programmable positions and four power functions in Hand and Foot Controls. Return-to-Home feature is designed for uninterrupted Patient-Practitioner Communication. Programming is a simple touch of a button and takes less than 2 seconds.

Patient-Centric Design

Double articulating headrest with height adjustability is designed for precise and comfortable head positioning for patients of all heights. Patient Weight Capacity of 650 lbs. accommodates a changing patient population, and ADA height of 18.5" enables seamless wheelchair transfers. The compact footprint makes this Procedure Chair the optimal choice for every setting.



DOUBLE-ARTICULATING HEADREST



OneTouch Locks and Unlocks the base. Glide your procedure chair around with ease.

The OneTouch WheelBase® is a game changer!

Integrated central-locking wheelbase enables staff to easily move the chair **for cleaning** and lock it solidly into place for exams and procedures.

ADA height
650 lb capacity
New compact design
OneTouch WheelBase®
Premium comfort cushion



Shown with optional Armrests (047)

“Vulnerable patient populations rely on frequent and intensive use of outpatient care to maintain or improve their health... It is critical that all of this care be provided under conditions that minimize or eliminate risks of healthcare-associated infections (HAI). **All healthcare settings, regardless of the level of care provided, must make infection prevention a priority.**”

Guide to Infection Prevention For Outpatient Settings: Minimum Expectations for Safe Care, CDC, National Center for Emerging and Zoonotic Diseases. Nov. 2015



The **proglide4010 Head Centric Procedure Chair** is designed to optimize the patient and caregiver experience. With ultimate accessibility to the head and neck area, the proglide 4010 by UMF Medical is setting a new benchmark for infection control, smart efficiencies, optimal accessibility, and lasting value. By asking what's going to improve your patients' and practitioners' experience, UMF Medical has arrived at a procedure room solution that is the perfect blend of style and function.

THIS IS THE PROCEDURE CHAIR YOU'VE BEEN WAITING FOR!

UMFmedical

1316 Eisenhower Boulevard
Johnstown, PA 15904-3393
Phone: 814-266-8726
Toll Free: 800-638-5322
Fax: 814-266-1870
www.UMFmedical.com

