



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Midmark 626 Power Exam Table Refurbished



Stock Photo – Options May Vary

The Midmark 626 Powered Exam Table offers exceptional patient accessibility, with a low height that allows easy entry and exit. It features powered back movement, adjustable pelvic tilt, and a heated drawer for enhanced comfort. The table is equipped with a hospital-grade power outlet, paper roll holder, and easy-to-use hand and foot controls for precise patient positioning.

Designed with patient comfort in mind, the table accommodates up to 650 lbs, making it ideal for bariatric patients. Additional features include a barrier-free, low chair height and a Clean Assist Roller System for effortless cleaning.

The Midmark powered exam table comes fully refurbished with new upholstery in your choice of colors, hand control, and a comprehensive warranty.

- **Easy Positioning:** Improves blood pressure accuracy.
- **Powered Back Movement:** Supports optimal back positioning.
- **Barrier-Free Height:** Fits most patients, feet flat on floor.
- **Clean Assist System:** Easy, safe chair movement for cleaning.
- **Hand & Foot Control:** Quick positioning at your fingertips.
- **Pelvic Tilt:** Ideal for lower body exams.
- **Active Sensing:** Pauses movement if disrupted.

Specifications:

- **Patient Weight Capacity:** 650 lb (295 kg)
- **Height:** Uncompressed: 15.5" (Low) / 37" (High) [39 cm / 94 cm] Compressed: 14.6" (Low) [36 cm]
- **Length (with Footrest Extended):** 76" (193 cm)
- **Pelvic Tilt:** 5.2°
- **Heated Drawer:** 24.6-Watt, adjustable from 92°F to 112°F.
- **Paper Roll Holder:** Fits 24" x 3.56" roll. **Electrical Receptacles:** Hospital-grade on both left and right sides, 4 A electrical draw
- **Table Weight:** 435 lb (197 kg)

***Brochure and images are a general overview of the product. Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2025 Venture Medical ReQuip, Inc.