

# Hausted® Mammography/Biopsy Chair

The Perfect Chair for Radiology  
and Mammography



## Stretcherology: MBC000ST

### Functional Benefits

- Designed specifically for Mammography
- Its unique layout allows for a full upright position to maneuver between the x-ray machine and the wall
- Self storing side rail allows the chair to maneuver into tight spaces
- Smallest footprint of the Hausted® chairs

### Versatility

- Can be used for pre-induction, transport, treatment and recovery procedures
- Trendelenburg and the emergency release lever equip the chair for immediate response to hemodynamic instability
- The pneumatic assisted back rest can adjust for supine to 90° for procedural positioning

### Patient Comfort

- High density, non-conductive vinyl covered foam pads support the patient comfortably
- Leg section and height of the footrest are adjustable, maximizing lower body comfort and stability



### Product Specifications

MODEL	MBC-000- ST
Seat Cushion Width	22 in [55 cm]
Overall Length (Supine)	75 in [190 cm]
Overall Length (Upright)	42.5 in [108 cm]
Height Adjustment	23.5 - 31.5 in [60 - 80 cm]
Wheel Base	25 x 22 in [63 x 55 cm]
Restraint Strap	Included
Pneumatic-Assist Backrest	0 - 90°
Trendelenburg	10°
Maximum Weight Capacity	325 lbs [147.5 kg] Evenly Distributed



770.368.4700

[www.hausted.com](http://www.hausted.com)

GF1400048RevA14



**GRAHAM-FIELD®**  
**Proudly made in the USA**

© 2014 GF Health Products, Inc